

Systech Illinois offer a range of adaptors to be used with the 7000 and 8000 range, 8200 and 8701 oxygen permeation analysers. Custom made adaptors are also available.



## Adaptor Glue Seal with Bracket

Glue any bottle or package with a neck opening of more than 15mm to the platen.



## Adaptor Film Test

This adaptor is used to perform standard film tests using 50cm<sup>2</sup> active area.



## Adaptor 1/8" NPTF

Using 1/8" NPT couplings and pipes the user can connect to any package in the traditional manner.



## Adaptor Complete Package

The simplest method to test a closed package. Drill a 1/4" hole in the package and glue the tube. Quick connect the assembly and start to test.



## Accessories & Spares

	Calibration/ Reference Film	Connect Tube	Coupling	Ferrule	Film Sealant	Male Connector	Metal Gasket	Nut	Tube	Stainless steel pipe
<p><b>Adaptor Glue Seal with Bracket</b></p> <p>Use 1/16" tube with fittings for flow inlet to package. Metal gasket needed for new connection to platen.</p> <p>7000 &amp; 8000 Part Number: 081 502 8200 &amp; 8701 Part Number: 870 051</p>				101 755			101 776 (10)	101 753	101 754 (1mtr)	
<p><b>Adaptor 1/8" NPTF</b></p> <p>Use 1/8" male connector and pipe to glue to any package. Metal gasket needed for new connection to platen.</p> <p>7000 &amp; 8000 Part Number: 081 503 8200 &amp; 8701 Part Number: 870 054</p>			101 757 (1/8" NPTF)			101 065 (1/8" NPT)	101 776 (10)			101 711 (1/8" 10mtr)
<p><b>Adaptor Complete Package</b></p> <p>Use 1/16" tube with fittings for flow inlet to package. Select a connect tube. Metal gasket needed for new connection to platen and for connect tube.</p> <p>7000 &amp; 8000 Part Number: 081 504 8200 &amp; 8701 Part Number: 870 062</p>		870 063 (150mm)  870 064 (250mm)		101 755			101 776 (10) (to platen)  870 059 (to connect tube)	101 753	101 754 (1mtr)	
<p><b>Adaptor Film Test</b></p> <p>Use with either Film sealant of high quality (good barriers) or low quality (poor barriers). Or with Calibration/Reference film, sample size A4, sufficient for two tests. Metal gasket needed for new connection to platen.</p> <p>8200 &amp; 8701 Part Number: 870 055</p>	080 080 (A4 size)				080 055 (high quality)  085 055 (low quality)		101 776 (10)			